

# Application of force and inertial sensors to monitor gait on legacy walkers <sup>†</sup>

Vítor Viegas<sup>1,2,3</sup>, J. M. Dias Pereira<sup>1,2</sup>, Octavian Postolache<sup>1,4</sup>, Pedro Silva Girão<sup>1,5</sup>

- <sup>1</sup> Instituto de Telecomunicações, 1049-001 Lisboa, Portugal
- <sup>2</sup> ESTSetúbal, <sup>3</sup>CDP2T, Instituto Politécnico de Setúbal, 2910-761 Setúbal, Portugal
- <sup>4</sup> ISCTE Instituto Universitário de Lisboa, 1600-077 Lisboa, Portuga
- <sup>5</sup> Instituto Superior Técnico, Universidade de Lisboa, 1049-001 Lisboa, Portugal
- <sup>†</sup> This paper is an extended version of our conference paper: Vítor Viegas, J. M. Dias Pereira, Octavian Postolache, Pedro Silva Girão,
- "Application of Force and Inertial Sensors to Monitor the Usage of Walker Assistive Devices", 22<sup>nd</sup> IMEKO TC4 Internation Symposium (2017)

#### **ABSTRACT**

Walker assistive devices play an important role in extending the autonomy of elderly people and in recovering the mobility of people affected by locomotion disabilities. The next generation of walkers will include embedded sensors and data processing capabilities that will allow the extraction of objective metrics, about gait and body posture, to assist the work of the physiotherapist and to enable the self-control of walker usage. The paper presents the *Andante*, our last proposal of a smart walker intended to monitor and analyze gait in real time. The system makes use of e-textile electrodes to sense the heart rate of the user, load cells to measure the forces applied on the walker legs, and an inertial measurement unit to sense motion and orientation. These signals are sampled locally and transferred over a Bluetooth link to a remote host where they are processed in real time. Data processing includes the detection, classification and characterization of the steps. A rich set of parameters is presented for each step, including estimates of balance and motor coordination, step length and azimuth, and lift of the walker frame. This information can be used by the physiotherapist to assess objectively the physical condition of the user and tune the rehabilitation therapy if needed. The proposed solution can be easily integrated into any commercial walker without any loss of functionality.

Section: RESEARCH PAPER

Keywords: force measurement; load cell; IMU; walker assistive device; gait analysis.

Citation: Vítor Viegas, J. M. Dias Pereira, Octavian Postolache, Pedro Silva Girão, "Gait monitoring on instrumented walkers: a proposal based on force and inertial sensors", Acta IMEKO, vol. 3, no. 1, article 1, January 2014, identifier: IMEKO-ACTA-03 (2014)-01-01

Editor: Paolo Carbone, University of Perugia, Italy

Receivedmonthday, year; In final formmonth day, year; PublishedJanuary2014

**Copyright:** © 2014IMEKO. This is an open-access article distributed under the terms of the Creative Commons Attribution 3.0 License, which permits unrestricted use, distribution, and reproduction in any medium, provided the original author and source are credited

Funding: This work was supported by Fundação para a Ciência e Tecnologia, project PTDC/DTP-DES/6776/2014.

Corresponding author: Vítor Viegas, e-mail: vitor.viegas@estsetubal.ips.pt

## 1. INTRODUCTION

Assistive walking devices play an important role in extending the autonomy of elderly people and in preventing accidental falls. This is especially important in modern societies where population is aging fast: in North America, 25% of the population will be over 65 years old by the year 2025 [1-2]; and in the UE, by the year 2060, the life expectance for women and men will be around 89 and 84.5 years, respectively [3].

Assistive walking devices are also key elements in recovering the mobility of people affected by locomotion disabilities (due to amputation, injuries of spine, muscular dystrophies or other causes). The use of walking aids is in

general positive because it improves the confidence levels during long-term rehabilitation therapies. A study published by the Vogt et al. [4] showed that the use of rollators (walkers with wheels) does not interfere with the rehabilitation outcome and, in some cases, may be very useful to decrease the rehabilitation periods.

It is consensual that the proper use of mobility aiding devices contributes to the sustainability of heath care systems [5-6], but the very same devices, if badly used, can be harmful and cause serious injuries [7-8]. The safety margin could be improved if users were educated on how to use the equipment properly, and if intelligence was incorporated on the devices to detect and avoid risk situations.

In many rehabilitation processes, the evaluation of user's balance and motor coordination is usually based on subjective data acquired through human observations. The decisions of the physiotherapist could be improved if they were based on objective metrics obtained from sensed data.

To address these problems, several instrumented walkers have been proposed to monitor and analyze gait [9-10]. Load measurements are usually done by placing force sensors on the walker hand grips or on the walker legs [11-15]. Kinematic parameters related with human gait are often measured with accelerometers or inertial motion units (IMU) attached to the walker frame. In this case, long-term drifts and persistent DC components must be filtered out to not affect the time integration of the signals [16].

Solutions based on RADAR measurements are also referred in the literature [17-18], and they are successful in which concerns the identification of gait patterns, but not so efficient to obtain accurate measurements about walking distance measurements.

Ultrasounds can also be a solution to obtain distance measurements [19], but they are not directional and are affected by environmental factors such as temperature and humidity.

Several works [20-21] also refer the usage of image processing techniques to monitor gait. The main advantages of these methods are related with their non-invasiveness, but there are strong limitations related with camera positioning and it is not possible, at all, to acquire data concerning force intensities.

Our team has been working on smart walkers over the last years [22-25]. We have tried a vast set of sensors, including force sensing resistors and load cells to measure force, and LIDAR and microwave Doppler RADAR to measure distance and velocity. Our work has two main purposes: first, monitor (and eventually correct) the usage of assistive gait devices; and second, extract objective gait metrics to assist physiotherapists in their work.

In this paper we present the Andante, our last proposal of a smart walker and the evolution of the prototype presented in [25]. It includes e-textile electrodes to capture the heart rate of the user, load cells to measure the forces applied on the walker legs, an IMU to sense motion and orientation, and software to bind all things together. Care was taken to preserve the native functionality of the device and to reduce the upgrade costs.

The paper is organized as follows: section II defines the geometry and metrics used along the paper; section III describes the implemented measurement system; section IV presents experimental results; and section V draws conclusions.

## 2. PRELIMINARY DEFINITIONS

A Cartesian coordinate system is established to represent the walker frame, and two indicators are defined to characterize user's gait.

# 2.1. Balance

The coordinate system illustrated in figure 1 will be used to study the forces applied on the walker frame. The walker legs are numbered from 1 to 4 (as quadrants) and the y-axis points to the forward direction. According to this arrangement, the centre of forces (COF) is given by:

$$COF_{x} = \frac{W_{12}(F_{1} - F_{2}) + W_{43}(F_{4} - F_{3})}{2(F_{1} + F_{2} + F_{3} + F_{4})} \tag{1}$$

$$COF_{\chi} = \frac{W_{12}(F_1 - F_2) + W_{43}(F_4 - F_3)}{2(F_1 + F_2 + F_3 + F_4)}$$

$$COF_{y} = \frac{L(F_1 - F_4) + L(F_2 - F_3)}{2(F_1 + F_2 + F_3 + F_4)}$$
(2)

where Fk represents the force applied on the kth leg, W12 and W43 represent the distances between the front and rear legs respectively, which corresponds roughly to the mean walker width (W), and L represents the distance between the front and rear legs, which corresponds to the walker length. In the present case, we are using a walker from Ayudas Dinámicas, model AD230 [26], with dimensions W12 = 51 cm, W43 = 53cm, L = 45 cm and weight = 2.5 kg.


The relative deviation of the COF in relation to the geometrical centre of the walker polygon is a measure of unbalance. Conversely, the complement to one is an indicator of balance (B):

$$B(\%) = 100 \times \left(1 - \frac{\sqrt{(COF_x)^2 + (COF_y)^2}}{\sqrt{(W/2)^2 + (L/2)^2}} \times \alpha\right)$$
(3)

Alpha (a) is a weighting factor that gives significance to B only when the user is loading the walker. It takes values from 0 to 1 depending if the walker is resting or is fully loaded with the

$$\alpha = \frac{F_1 + F_2 + F_3 + F_4}{F_{II}} = \frac{F_T}{F_{II}} \tag{4}$$

where F<sub>T</sub> represents the total force applied on the walker legs and F<sub>U</sub> represents the user weight.






Figure 1. Coordinate system of the walker frame: a) instrumented walker; b) Cartesian plane, dimensions and coordinates.




Figure 2. Gait phases (GP) during a walker step (the polygon delimits the support area). GPO: the walker is resting on the floor; GP1: the walker is flying; GP2: the walker is on the floor waiting for the injured foot to move forward; GP3: the injured foot is moving forward; GP4: waiting for the healthy foot to move forward; GP5: the healthy foot is moving forward.

## 2.2. Motor coordination

Motor coordination has to do with the right sequence of movements needed to complete a step. Considerer the state machine illustrated in figure 2, where the states represent the gait phases and the continuous arrows represent normal transitions between states [27-28]. If the user completes a step passing through all the gait phases, the step is classified as *good*; otherwise is classified as *bad*. The number of *good* steps (G) in the last N is an indicator of motor coordination (MC):

$$MC(\%) = 100 \times \frac{G}{N} \tag{5}$$

## 3. MEASUREMENT SYSTEM

The measurement system includes a wireless data acquisition module, an IMU, a heart rate sensor, four load cells and software to bind all things together.

# 3.1. Data acquisition

Data acquisition is performed by a Shimmer3 module from Shimmer Sensing [29]. The module includes an embedded IMU, four external analog inputs (AI1-AI4) and a Bluetooth link. Wiring was passed through the tubular structure of the walker to connect the sensors to the analog inputs. Samples from all sensors (IMU and external) are acquired at a rate of 51.2 S/s and sent wirelessly to a remote host. The sampling rate is well beyond the expected Nyquist frequency for physical movements.

# 3.2. Inertial sensors

The embedded IMU is used to sense motion and orientation. It has nine degrees of freedom (9DoF) given by a 3-axes accelerometer (model KXRB5-2042 from Kionix), a 3-axes gyroscope (model MPU-9150 from Invensense) and a 3-axes magnetometer (model LSM303DLHC from STMicroelectronics).

The IMU is placed horizontally with its coordinate system aligned to that of the walker (shown in figure 1). This way, the acceleration across the y axis gives motion in the forward-backward direction, and the acceleration across the z axis, subtracted by  $1 g = 9.81 \text{ ms}^{-2}$ , gives motion in the vertical.

The Shimmer3 module is also able to know its orientation by fusing the data provided by the IMU in a so-called "gradient descent algorithm" [30]. The algorithm gives orientation in absolute ENU coordinates (local East-North-Up) with a delay of few tenths of a second.

The IMU must be calibrated in advance to cancel offset errors and tune gain matrices. The calibration procedure [31] is executed once, in situ, and the results are saved into a non-volatile memory.

#### 3.3. Heart rate sensor

The heart rate is sensed by a pair of dry e-textile electrodes placed over the handgrips of the walker, tied with a tiny Velcro strip. They are removable, washable and seem like ordinary cloth, thus not causing additional stress to the user.

The ECG signal is conditioned by an analog front-end based on the AD8232. The signal passes through a two-pole band-pass filter, with cut-off frequencies at 7 Hz and 24 Hz, followed by an amplifier with gain 1000. The narrow-band filter eliminates motion artefacts and line noise but distorts the ECG waveform significantly, making it suitable only for heart-rate detection. The output signal is connected to the first analog input (AI1) of the data acquisition module.

# 3.4. Force sensors

The forces applied on the walker legs are sensed by four load cells. Each cell is attached to the extremity of a leg using a dedicated plastic adapter grown on a low-cost 3D printer, as shown in figure 3a. The cell contains four strain gauges fixed to a small bending beam (55.3 x 12.7 x 12.7 mm) that supports 20 kgf. Other characteristics include: nominal output sensitivity = 1 mV/V, non-linearity = 0.05% FS (full-scale) and hysteresis = 0.05% FS. Bending beam cells were preferred to inline/axial cells because they are much cheaper.

Some engineering had to be done to fit four load cells into the three analog inputs available:

- Consider the auxiliary axes p and q, passing over the load cells 1 & 3, and 2 & 4, respectively, as shown in figure 3b. The pair (p,q) is not perfectly orthogonal because the mean walker width (W) is not equal to the walker length (L). Nevertheless, the pair is *almost* orthogonal, with a rotation angle *near* to +45 degrees relatively to (x,y).
- Each load cell is supplied at 3 V and conditioned by an instrumentation amplifier with gain 100. The voltages are combined in pairs, subtracted from each other (1-3 and 2-4), amplified by 4 and biased around 1.5 V, as shown in figure 3c. The result is a voltage that varies linearly between 0.3 V and 2.7 V for forces applied along the axes p and q, between -D/2 and +D/2, where D = √L² + W² is the walker diagonal. The COF coordinates over the auxiliary axes are given by:

$$COF_p = \left(\frac{D}{2.4}\right) \times (V_{AI2} - 1.5) \tag{6}$$


$$COF_q = \left(\frac{D}{2.4}\right) \times (V_{AI3} - 1.5) \tag{7}$$


where  $V_{\Lambda IX}$  represents the voltage acquired at the analog input 2 or 3, in volt.


• The conversion from (p,q) to (x,y) coordinates is done by applying a rotation matrix with  $\theta = -45$  degrees:

$$\begin{bmatrix} COF_x \\ COF_y \end{bmatrix} = \begin{bmatrix} \cos \theta & \sin \theta \\ -\sin \theta & \cos \theta \end{bmatrix} \begin{bmatrix} COF_p \\ COF_q \end{bmatrix}$$
(8)

 Additionally, the signals from all load cells are summed together and injected in the analog input AI4 to get the total load applied on the walker. The resulting voltage varies linearly between 0 and 1.2 V for a total load between 0 and 80 kgf, thus giving:







LCX: Load Cell X; Blocks based on: a) AD8232; b) AD623; c) MCP6272.

Figure 3. Load cells attached to the walker legs: a) photo; b) auxiliary axes; c) signal conditioning.

$$F_T = \left(\frac{80}{12}\right) \times V_{AI4} \tag{9}$$

Playing with the walker symmetry, together with some basic geometry, we were able to fit everything in just three analog inputs.

# 3.5. Software

The samples arriving to the host are handled by the *Andante* application, which acts as the interface to the physiotherapist. The application was developed in Visual Studio 2012 and makes use of the Shimmer C# API [32] to interact with the Shimmer3 module.

The Andante application records its state on a relational database designed in Microsoft Access. The database makes possible to see how things evolve over time (e.g. query the performance of a given user during several therapy sessions) and to trace actions (e.g. query which therapist assisted a given user on a given session). The database is composed by five related tables as shown in figure 4:

- Physiotherapists: List of physiotherapists registered in the application. Each physiotherapist has a unique identifier (ID) associated to its name.
- Users: List of users registered in the application. Each user is described by ID, name, age, weight and injury (left or right limb).
- Walkers: List of four-leg instrumented walkers registered in the application. Each walker is characterized by ID, brand, model, serial number, width between front legs (W12), width between rear legs (W43), distance between front and rear legs (L) and offset errors.
- Sessions: List of physiotherapy sessions supervised by the application. Each session is characterized by ID, physiotherapist in charge, user under treatment, walker in use, timestamp and location of session opening.
- Steps: List of steps detected during a given physiotherapy session. Each step is characterized by ID, gait metrics computed by the measurement system and free comments added by the physiotherapist.

The tables of physiotherapists, users and walkers must be edited in advance, out of the *Andante* application, using Microsoft Access© tools for instance. The other tables are updated by the *Andante* application in runtime.




Figure 4. Main window of the *Andante* application.

The main window of the *Andante* application has five buttons that give access to the following options (see figure 5):

- Login: Opens a dialog where the physiotherapist enters the data needed to create a new physiotherapy session (physiotherapist in charge, user under treatment and walker in use).
- Connect: Establishes a Bluetooth connection with the instrumented walker.
- Calibrate Sensors: Opens a dialog to extract the offset errors of the measurement system. With the walker resting on the floor, with no external forces applied, the physiotherapist is instructed to register the COF coordinates and the weight of the walker frame. These residues are stored in the database, in the record corresponding to the walker, columns COFx0, COFy0 and F0. Later, during data streaming, these residues are subtracted from equations (1), (2) and from the numerator of equation (4).
- Stream Data: Opens a child window where the physiotherapist can monitor the usage of the walker (see section 3.5.1).
- Disconnect: Closes the Bluetooth connection with the instrumented walker.




Figure 5. Main window of the Andante application.

# 3.5.1. Stream data child window

The Stream Data child window opens a high-speed communication link to transfer data from the instrumented walker to the host (the link is unidirectional). The window launches two threads: a high-priority thread that stores the fast-arriving samples into a buffer; and a low-priority thread that runs every 200 ms to process the buffered data and update the user interface. The buffer contains a variable number of samples, from 10 to a few more, depending on the computer load. An update rate of 5 frames/s is enough to give the application good responsiveness and fluidity.

The buffered samples are processed sequentially, by the same order as they arrived. The processing includes the following main tasks:

- Computation of the heart rate from the signal connected to AI1. A simple Schmitt-trigger algorithm, with fixed thresholds at 0.5 V and 1 V, is executed to detect the heart beats. The time between beats is used to compute the heart rate.
- 2. Computation of COFx and COFy from the signals connected to AI2 and AI3. This is done by solving equations (6), (7) and (8), in this order.

- 3. Computation of the balance indicator (B) considering the signal connected to AI4. This is done by solving equations (9), (4) and (3), in this order.
- 4. Detection of *good* steps by running the state machine presented in figure 6.
- 5. Computation of the motor coordination indicator (MC) by solving equation (5).

The state machine illustrated in figure 6 is a key element of the instrumented walker because it classifies steps as *good* or *bad*. States 0 to 5 represent the gait phases (as in figure 2) and state 10 represents a waiting state to rearm the machine. Transitions between states ( $T_{from,to}$ ) are triggered by sensor measurements as follows:

- T<sub>0,1</sub> is enabled if F<sub>T</sub> < WEIGHT\_TH. This transition occurs when the total force measured by the load cells falls below a given threshold (below the walker weight). In other words, it occurs at the beginning of a new step when the walker is lifted in the air.
- $T_{1,2}$  is enabled if  $F_T$  > WEIGHT\_TH. This transition occurs when the walker touches down the floor again.
- T<sub>2,3</sub> is enabled if COFx > RIGHT\_TH or COFx < LEFT\_TH depending on the injury (left or right limb, respectively). This transaction occurs when the user moves forward the injured limb and applies force on the opposite side. The injured limb shall be the first to move forward. The type of disability must be defined in advance to determine the threshold value.
- T<sub>3,4</sub> is enabled if COFx < RIGHT\_TH or COFx > LEFT\_TH depending on the injury (left or right limb, respectively). This transaction occurs when the user




Figure 6. State machine used to classify steps as good or bad.

alleviates the force previously applied on the walker making COFx return to the origin.

- T<sub>4,5</sub> is enabled if COFx < LEFT\_TH or COFx > RIGHT\_TH depending on the injury (left or right limb, respectively). This transaction occurs when the user moves forward the healthy limb and applies force on the opposite side. The healthy limb shall be the last to move forward.
- T<sub>5,0</sub> is enabled if COFx > LEFT\_TH or COFx < RIGHT\_TH depending on the injury (left or right limb, respectively). This transaction occurs when COFx returns to the origin.
- T<sub>2,0</sub> and T<sub>4,0</sub> are enabled if COFx moves to the wrong side, i.e. if the user moves forward the wrong limb. The machine returns to state 0 and waits for a new step.
- T<sub>2,10</sub> and T<sub>4,10</sub> are enabled if F<sub>T</sub> < WEIGHT\_TH. This transition occurs if the user lifts the walker before finishing the step. In other words, it occurs if the user aborts the step.
- T<sub>10,0</sub> is enabled if F<sub>T</sub> > WEIGHT\_TH. This transition occurs when the walker touches down the floor again. The machine returns to state 0 and waits for a new step.

A step is marked as *good* if the machine passes through all the states from 0 to 5. A step is marked as *bad* if the user aborts it or moves forward the wrong limb.

The threshold WEIGHT\_TH is equal to half of the weight of the walker frame. The thresholds RIGHT\_TH and LEFT\_TH are dynamic and computed as follows:

(10) 
$$\begin{cases} LEFT_{-}TH(0) = -\frac{W}{6} \\ LEFT_{-}TH(k) = \frac{minCOF_{x}(k-1)}{2} \\ If (LEFT_{-}TH(k) > -50 \text{ mm}) LEFT_{-}TH(k) = -50 \text{ mm} \end{cases}$$

(11) 
$$\begin{cases} RIGHT_{-}TH(0) = +\frac{W}{6} \\ RIGHT_{-}TH(k) = \frac{maxCOF_{x}(k-1)}{2} \\ If (RIGHT_{-}TH(k) < +50 \text{ mm}) RIGHT_{-}TH(k) = +50 \text{ mm} \end{cases}$$

where k represents the kth *good* step, and min/maxCOFx(k-1) represents the minimum/maximum value reached by COFx during the previous (k-1) *good* step.

The processed data is shown in the graphical interface shown in figure 7, which contains four panels:

- Panel 1: Shows the instantaneous value of the heart rate.
   This is useful to estimate the effort and stress felt by the user.
- Panel 2: The COF is computed, and the result is presented as a cross moving over a XY graph. When the user loads his left side the cross moves toward negative values of X; otherwise, the cross moves toward positive values of X. The same applies for the front/back direction over the Y axis, much like a joystick. A vertical slider gives the instantaneous value of balance (B).
- Panel 3: A set of six LEDs are turned on sequentially while the state machine moves forward from state 0 to state 5. If the user passes through all the gait phases successfully, all the LEDs light up and the step is marked as *good*. If the user violates any gait phase, the machine is reset to state 0 and the step is marked as *bad*. The number of *good* and *bad* steps is registered. A vertical slider gives the instantaneous value of the motor coordination (MC).
- Panel 4: Presents a list of all steps detected. For each step the following information is provided:
  - Step#: Index of the step.
  - Timestamp: Date and time of the step.
  - Heart rate: Last value of heart rate measured during the step.




Figure 7. Graphical interface of the Stream Data child window.

- Balance: Minimum value of balance (B) measured during the step.
- Motor coordination: Value of motor coordination (MC) considering the current step.
- Quality: Classification of the step (good or bad).
- Lift: Lift of the walker frame while moving forward, in millimetre. The lift is estimated by double integrating the acceleration across the z axis, subtracted by 1 g = 9.81 ms<sup>-2</sup>. Prior to integration, the acceleration samples, collected only during state 1, are smoothed to remove high-frequency noise, and are subtracted from the baseline to remove drift. These two actions are equivalent to apply a bandpass filter.
- Step length: Length of the step, in millimetre. The
  distance is estimated by double integrating the
  acceleration across the y axis. Prior to integration, the
  acceleration samples, collected only during state 1, are
  bandpass filtered as described above.
- Azimuth: Azimuth of the walker frame, in degrees. The values {0°, 90°, 180°, 270°} correspond to {north (magnetic), east, south, west}, respectively.
- Failure: Details de quality of the step according to the following enumeration: none, step aborted, injured foot failed to move forward or healthy foot failed to move forward.

The physiotherapist can also add free comments to each step. The list can be saved to the database (into the Steps table) by clicking on a button.

## 4. EXPERIMENTAL RESULTS

The experimental tests were made using a standard four-leg walker, model AD230 from Ayudas Dinámicas [26], with dimensions W12 = 51 cm, W43 = 53 cm, L = 45 cm and weight = 2.5 kg.

The first test was intended to evaluate the accuracy of COF measurements. The four walker legs were placed over four digital scales, model Spider 1S-60S from Mettler Toledo, and load was applied on the walker frame to move the COF to the 2<sup>nd</sup> quadrant. The COF coordinates were computed knowing the geometry of the walker and the weights indicated by the scales (with a resolution of 20 gf). The COF coordinates measured by the *Andante* application were compared against the COF coordinates computed from the scales (the scales serve as reference).

The experiment was repeated rearranging the load to move the COF to the  $4^{th}$  quadrant. Table 1 shows the COF coordinates measured by both systems (*Andante* and scales) in both arrangements ( $2^{nd}$  and  $4^{th}$  quadrants). The distances to the origin (r) are also shown, as well as the corresponding relative errors ( $\eta$ ). The values obtained, below 10%, are acceptable considering that the load cells are low cost and their alignment is not perfectly orthogonal (as explained in section 3.4).

Table 1. Comparison between the *Andante* and a reference system regarding the computation of COF coordinates.

Experiment	COF coordinates		Distance to origin r		
	Andante	Scales	Andante	Scales	η
2 <sup>nd</sup> quadrant	(-122; 78) mm	(-126; 96) mm	145 mm	158 mm	8.2%
4 <sup>th</sup> quadrant	(180; -44) mm	(194; -45) mm	185 mm	199 mm	7.0%

The second test was intended to see how the COF swings to the right and to the left during a *good* step. Figure 8 shows precisely that, together with the value of balance (B). Both are related, with the balance falling whenever the COF moves away from the origin. The curves are explained as follows:

- State 0 was omitted because it is a waiting state.
- State 1 begins when the walker is lifted in the air and ends when the walker touches the floor again.
- State 2 is a waiting state that ends when the user moves forward the injured limb (making COFx cross LEFT TH).
- During state 3, the user loads his left side to compensate the lack of support on the injured limb. Then COFx returns to the origin and the machine moves to the next state.
- State 4 is another waiting state that ends when the user moves forward the healthy limb (making COFx cross RIGHT\_TH).
- During state 5, the user transfers part of his weight to the right. The load applied during state 5 is lower than that applied on state 3 because the healthy limb moves more easily than the injured limb. The step ends when both feet stay side-by-side, making COFx return to the origin and restart the state machine.





Figure 8. Synchronized plots of COFx and balance (B) during a *good* step. The numbered circles have correspondence with the states of figure 6.

The third test was made with the help of a user with impaired gait caused by an injury in the right lower limb. The user was instructed to follow an imaginary straight line, first in an early stage of the therapy (figure 9a), and some weeks later in a later stage of the therapy (figure 9b). The XY graphs cover 12 steps, with each step being represented by a line segment. The length and angle of the segment represent the travelled distance and azimuth of the walker. The user moves from the origin (0,0) to the periphery.

Table 2 summarizes the metrics obtained for both experiments. The averages and standard deviations cover the 12 steps, except for the motor coordination that considers only the final value.

The user followed a trajectory almost linear in both experiments, with mean squared errors equal to 9.3 cm and 2.2 cm, in the early and final stages of the therapy, respectively. In fact, all metrics show that the physical condition of the user has improved: balance and motor coordination increased, the total travelled distance also increased, steps became longer and more regular, and the walker began to be handled more easily (as indicated by higher lifts).






Figure 9. User following a linear trajectory: a) in the early stage of the therapy; b) in a later stage of the therapy.

Table 2. Gait metrics acquired while the user was following an imaginary straight line.

		Experiment	
Quantity	Sub-quantity	Early stage	Final stage
Balance B (%)	average	49	78
MC (%)	final	60	90
Malkor lift (om)	average	1.3	6.5
Walker lift (cm)	std	1.2	1.5
Step length (cm)	average	21.2	33.6
step length (tm)	std	9.4	8.9
Total travelled distance (m)		2.54	4.03

# 5. CONCLUSIONS

The experimental results show that the *Andante* is able to detect and classify steps on a standard four-leg walker. Each step is characterized by a rich set of parameters, including estimates of balance and motor coordination, step length and azimuth, and lift of the walker frame. This record of information can be used by a physiotherapist to assess objectively the physical condition of the user. If these records are saved during multiple sessions, the physiotherapist can look at them, understand how the rehabilitation therapy is progressing over time, and take corrective actions if needed.

## **ACKNOWLEDGEMENT**

The work was supported by Instituto de Telecomunicações and Fundação para a Ciência e Tecnologia, projects PTDC/DTPDES/1661/2012 and PTDC/DTP-DES/6776/2014.

## **REFERENCES**

- [1] Carrie A. Werner. "The older population: 2010". 2010 Census Briefs, U.S. Census Bureau. <a href="https://www.census.gov/prod/cen2010/briefs/c2010br-09.pdf">https://www.census.gov/prod/cen2010/briefs/c2010br-09.pdf</a> (February 2018)
- [2] "Older Americans 2008: key indicators of well-being". Federal Interagency Forum on Aging-Related Statistics. https://agingstats.gov/docs/PastReports/2008/OA2008.pdf (February 2018)
- [3] "The demographic future of Europe from challenge to opportunity". European Commission. <a href="http://europa.eu/legislation\_summaries/employment\_and\_social\_policy/situation\_in\_europe/c10160\_en.htm">http://europa.eu/legislation\_summaries/employment\_and\_social\_policy/situation\_in\_europe/c10160\_en.htm</a> (February 2018)
- [4] L. Vogt, K. Lucki, M. Bach, W. Banzer. "Rollator use and functional outcome of geriatric rehabilitation". Journal of Rehabilitation Research & Development, Issue 47, No. 2, pp. 151-156, 2010
- [5] Sara M. Bradley, Cameron R. Hernandez. "Geriatric assistive devices". American Family Physician Journal, Issue 15, No. 84(4), pp. 405-411, August 2011
- [6] S. R. Faruqui, T. Jaeblon. "Ambulatory assistive devices in orthopaedics: uses and modifications". Journal of the American Academy of Orthopaedic Surgeons, Issue 18, No. 1, pp. 41–50, 2010
- [7] H. Bateni, B. E. Maki. "Assistive devices for balance and mobility: benefits, demands, and adverse consequences". Archives of Physical Medicine and Rehabilitation, Issue 86, No. 1, pp. 134–145, 2005
- [8] J. A. Stevens, K. Thomas, L. Teh, Greenspan et al. "Unintentional fall injuries associated with walkers and canes in older adults treated in U.S. emergency departments". Journal of the American Geriatrics Society, Issue 57, No. 8, pp. 1464–1469, 2009
- [9] Joel A. DeLisa et al. "Gait analysis in the science of rehabilitation". Rehabilitation research and development service, Department of veterans affairs, 1998
- [10] Alvaro Muro-de-la-Herran, Begonya Garcia-Zapirain, Amaia Mendez-Zorrilla. "Gait analysis methods: an overview of wearable and non-wearable systems, highlighting clinical applications". Sensors Review, Vol. 14, Issue 2, pp. 3362-3394, 2014
- [11] R. A. Bachschmidt, G.F. Harris, G.G. Simoneau. "Walker-assisted gait in rehabilitation: a study of biomechanics and instrumentation". IEEE Transactions on Neural Systems and Rehabilitation Engineering, Vol. 9, Issue 1, pp. 96 105, 2002
- [12] Majd Alwan et al. "Basic walker-assisted gait characteristics derived from forces and moments exerted on the walker's handles: results on normal subjects". Medical Engineering and Physics, No. 29, pp. 380-389, 2007
- [13] A. D. C. Chan, J. R. Green. "Smart rollator prototype". IEEE International Workshop on Medical Measurements and Applications (MeMeA), Ottawa, Canada, May 2008
- [14] J. M. Dias Pereira, Octavian Postolache, Vítor Viegas. "A low-cost measurement system to extract kinematic parameters from walker devices". IEEE International Instrumentation and Measurement Technology Conference (I2MTC), Pisa, Italy, 11-14 May 2015
- [15] E. Sardini, M. Serpelloni, M. Lancini, "Wireless instrumented crutches for force and movement measurements for gait monitoring". IEEE Transactions on Instrumentation and Measurement, Vol. 64, No. 12, pp. 3369-3379, December 2015
- [16] I-Min Lee, E. J. Shiroma. "Using accelerometers to measure physical activity in large-scale epidemiological studies: issues and

- challenges". British Journal of Sports Medicine, Vol. 48, Issue 3, pp. 197-201, February 2014
- [17] O. Postolache, P. Silva Girão, G. Postolache, G. Mendes. "Cardio-respiratory and daily activity monitor based on FMCW Doppler radar embedded in a wheelchair". Proceedings of 33<sup>rd</sup> Annual International Conference of the IEEE Engineering in Medicine and Biology Society, pp. 1917-1920, Boston, MA, USA, September 2011
- [18] G. Postolache, P. Silva Girão, O. Postolache. "Applying smartphone apps to drive greater patient engagement and personalized physiotherapy". Proceedings of IEEE International Symposium on Medical Measurements and Applications (MeMeA), pp. 1-6, Lisbon, Portugal, June 2014
- [19] A. Ekimov, J. M. Sabatier. "Human motion analyses using footstep ultrasound and Doppler ultrasound". Journal of the Acoustic Society of America, Vol. 123, No. 6, pp. 149-154, June 2008
- [20] Xi Cai, Guang Han, Jinkuan Wang. "Assessment of gait parameters during unconstrained straight-line walking using single camera". Proceedings of the IEEE 7th International Conference on Awareness Science and Technology (iCAST), September 2015
- [21] Xi Cai, Guang Han, Xin Song, Jinkuan Wang. "Single-camerabased method for step length symmetry measurement in unconstrained elderly home monitoring". IEEE Transactions on Biomedical Engineering, Vol. 64, Issue 11, pp. 2618-2627, January 2017
- [22] O. Postolache, P. Girão, J. M. Dias Pereira, G. Postolache. "Smart walker for pervasive healthcare". Proceedings of International Conference on Sensing Technology (ICST), Vol. 1, pp. 1-5, Palmerston North, New Zealand, December 2011
- [23] Octavian Postolache, J. M. Dias Pereira, Vítor Viegas, Pedro Silva Girão. "Gait rehabilitation assessment based on microwave

- Doppler radar embedded in walkers". IEEE International Symposium on Medical Measurement and Applications (MeMeA), Torino, Italy, 7-9 May 2015
- [24] J. M. Dias Pereira, Vítor Viegas, Octavian Postolache, P. Silva Girão. "Combining distance and force measurements to monitor the usage of walker assistive devices". IEEE International Instrumentation and Measurement Technology Conference (I2MTC), Torino, Italy, 22-25 May 2017
- [25] Vítor Viegas, J. M. Dias Pereira, Octavian Postolache, Pedro Silva Girão. "Application of force and inertial sensors to monitor the usage of walker assistive devices". 22nd IMEKO TC4 International Symposium, Iasi, Romania, September 2017
- [26] "Manual de usuário AD230AD/ AD233AD/ AD232". Ayudas Mecánicas. <a href="https://www.ayudasdinamicas.com/caminador-de-incorporacion.pdf">https://www.ayudasdinamicas.com/caminador-de-incorporacion.pdf</a> (February 2018)
- [27] David A. Winter, "Biomechanics of normal and pathological gait: implications for understanding human locomotor control", Journal of Motor Behavior, Vol. 21, No. 4, pp. 337-355, December 1989
- [28] Marder E., Calabrese R. L. "Principles of rhythmic motor pattern generation", Physiological Reviews, Vol. 76, Issue 3, pp. 687– 717, July 1996
- [29] <a href="http://www.shimmersensing.com/products/shimmer3-imu-sensor">http://www.shimmersensing.com/products/shimmer3-imu-sensor</a> (February 2018)
- [30] Sebastian O. H. Madgwick, Andrew J.L. Harrison, Ravi Vaidyanathan. "Estimation of IMU and MARG orientation using a gradient descent algorithm". IEEE International Conference on Rehabilitation Robotics, Zurich, Switzerland, June 2011
- [31] "9DoF Calibration Application User Manual Rev. 2.5a", Shimmer, 2014
- [32] http://www.shimmersensing.com/products/shimmercapture (February 2018)